


Dane aktualne na dzień: 29-05-2024 21:54

Link do produktu: <http://www.novazym.sklep.pl/sample-concentrator-nd100-1-p-4016.html>



Sample Concentrator ND100-1

Opis produktu

Sample Concentrator

ND100-1



Description:

- ◆ ND100-1 sample concentrator is with microprocessor controlled and PID fuzzy controlled technique. Instead of traditional water bath device, the heat-conducting medium is high purity aluminum. The Test tubes are held in a block, and blowing a large number of samples at one time. Can be used in residue analysis, commodity inspection, food, environment, pharmacy and biological product and other industries.

Features:

- ◆ High precision and wide range of temp control.
- ◆ Blowing a large number of samples at one time.
- ◆ Blowing each sample independently.
- ◆ Amount of airflow can be adjusted.
- ◆ Easy operation; safe and reliable to use.
- ◆ 12-month warranty.

Parameters:

Project	Parameters
Temp. Range	RT.+5 -160
Temp. Uniformity	±0.5
Display accuracy	0.1
Heating time	≤15 min (25 TO 160)
Max gas usage	15L/min
Max gas pressure	0.02Mpa (Gas pin ≤16pcs) 0.05Mpa (Gas pin >16pcs)
Max power	300W
Size	260*220*450mm

Project	Parameters
Temp accuracy	±0.5 (@40)
Temp accuracy	±1 (@120)
Time Range	1-99h59min/∞
Max vertical travel	200mm
Needle length	150mm
Block qty	1
Voltage	220V 50/60Hz
Weight	7.0KG

Optional Blocks:

Type	Test tube diameter	Capacity
MD01	6mm	48
MD02	7mm	48
MD03	10mm	24
MD04	12mm	24
MD05	13mm	16
MD06	15mm	16
MD07	16mm	16
MD08	19mm	12

Type	Test tube diameter	Capacity
MD09	20mm	12
MD10	26mm	8
MD11	28mm	4
MD12	40mm	3
MD13	0.2ml (centrifuge tube)	48
MD14	0.5ml (centrifuge tube)	48
MD15	1.5ml (centrifuge tube)	24
MD16	2.0ml (centrifuge tube)	24

Novazym
life science & technology

Novazym Polska S.C. Wielkopolska Centre of Advanced Technologies
ul. Umultowska 89 C build. A/38, 61-614 Poznań