

Link do produktu: <http://www.novazym.sklep.pl/ultrabright-led-transilluminator-dual-wavelengths-630-740-nm-p-388.html>



UltraBright LED transilluminator (dual wavelengths 630 & 740 nm)

Dostępność

Na zamówienie

Numer katalogowy

LB-17N

Opis produktu

[UltraBright LED transilluminator \(dual wavelengths\) MaestroGen](#)

Specification

- Illumination type: LED
- Wavelength: 740 nm (NIR), 630 nm (Red)
- Dimension (mm): 340L*280W*80H
- Viewing surface (mm): 160D*200W
- LED lamps arrangement: Delta matrix
- LED lifetime: 50,000 hrs
- Weight: 3.6 kg
- Power: 24 V - 1.25 A

[Key features](#)

- Application (In-Cell westerns) as an excitation source
- Strong R/NIR LEDs for excitation
- Dual wavelengths embedded in the same transilluminator
- Delta LED arrangement
- High uniformity across viewing area
- US patent pending
- CE/FCC Certified