

Link do produktu: <http://www.novazym.sklep.pl/dry-bath-incubator-mk2000-2-p-295.html>



Dry Bath Incubator MK2000-2

Dostępność

Na zamówienie

Numer katalogowy

MK2000-2

Opis produktu

Dry Bath Incubator MK2000-2

The MK2000-2 Dry bath incubator is ideal for incubation and activation of cultures, enzyme reactions, blood urea nitrogen determinations, immunoassays, melting/boiling points, and a wide variety of other laboratory procedures. The economical, superb temperature control and compact dry bath incubators are designed for applications that require repeatable results and superior temperature stability.

Specification

Temperature Control Range	RT+5 °C ~ 105 °C
Stability @ 100 °C	±0.4 °C
Uniformity:	
within the block @ 37 °C	±0.2 °C
across similar blocks @ 37 °C	±0.3 °C
Heating time (25 °C to 100 °C)	≤15min
Display Accuracy / Display	0.1 °C / LED
Timing Range	1min ~ 99h59min
Block Quantity	2
Dimensions L x W x H (mm)	365 x 210 x 150
Weight (kg)	4.5
Standard Blocks	See "Blocks for MK2000-2 Dry Bath Incubator"

Key features

- Exceptional temperature uniformity and accuracy -- due to advanced temperature control combined with high quality, precision engineered block providing excellent thermal contact
- Wide range of interchangeable aluminium blocks
- Heat to 105 °C, optional external temperature probe
- Blocks available as accessories for all applications - tubes, vials and microplates
- Block extraction tool is supplied, allowing blocks to be removed easily
- 2-years warranty as standard