

Link do produktu: <http://www.novazym.sklep.pl/bacti-cinerator-sterilizer-hm-3000-p-374.html>



# Bacti-Cinerator Sterilizer HM-3000

Dostępność

**Na zamówienie**

Numer katalogowy

**HM-3000**

## Opis produktu

### Bacti-Cinerator Sterilizer HM-3000

The Bacti-Cinerator Sterilizer sterilizes microorganisms utilizing infrared heat produced by a ceramic core element. The ceramic element contains no asbestos and ensures maximum sterility without splatter across the work surface area. Complete sterilization occurs within 5-7 seconds at optimum sterilizing temperature of 1 500 °F (815.6 °C). The small footprint makes the Bacti-Cinerator Sterilizer an ideal instrument for anaerobic and aerobic chambers alike. Several user friendly accessories are also available for the Bacti-Cinerator sterilizer.

#### Specification

1. The center high-temperature: 815 ±50 °C
2. The Max-diameter of sterilizer goods:  $\varnothing$  14 mm
3. The length of sterilizer: 140 mm
4. Net weigh: 1.3 kg

#### Key features

1. Sterilizes needles, loops, and culture tube mouths in five to seven seconds at optimum sterilizing temperature of 815.6 °C (1 500 °F)
2. Ceramic funnel tube enclosed in stainless-steel perforated guard and stainless-steel support stand
3. For AC 220 V, 50 / 60 Hz