

Link do produktu: <http://www.novazym.sklep.pl/biological-indicator-incubator-minit-1-p-518.html>



Biological Indicator Incubator MiniT-1

Dostępność

Dostępny

Opis produktu

Biological Indicator Incubator Model? MiniT-1

The MiniT-1 Biological Indicator Incubator is designed for temperature incubation of the EO Biological Indicators. It provides easy verification of incubation temperature.

Features

- 1. LCD display, with 12V power for ease of portability
- Operating Temperature Range: From 5°C above ambient to 80°C
- Time Range: 1min ~ 99h59min
- 4Temp.Control Accuracy: $\leq \pm 0.5^\circ\text{C}$
- 5Display Accuracy: 0.1°C
- 6Heating Time: < 6min (from 20°C to 57°C)
- Temp.Control Accuracy @57°C: $\leq \pm 0.35^\circ\text{C}$
- Dimensions: 110(w) x 150(D) x 80(H) mm

Specification

- 1. The Incubator has many wells for placement of the biological indicator vials and the wells can be custom-made. User can set the incubating temperature.
- 2. Materials: Aluminium block
- 3. Well Diameter: According to custom-making

Type	Description	Remark
A-MiniT-1	15 x ϕ 10mm	
B-MiniT-1	12 x ϕ 12mm	
C-MiniT-1	6 x ϕ 15mm	
D-MiniT-1	6 x ϕ 16mm	
E-MiniT-1	18 x ϕ 8.8mm	