

Dane aktualne na dzień: 31-05-2024 20:39

Link do produktu: <http://www.novazym.sklep.pl/dry-bath-incubator-dh100-2-p-4018.html>



Dry Bath Incubator DH100-2

Opis produktu

Dry Bath Incubator

DH100-2



Description:

- ◆ DH100 series high temperature dry bath incubator controlled by microcomputer, use of high purity aluminum material as a heat conduction medium instead of the conditional water bath device. which is easy to use, wide temperature range, high precision and etc. Can be widely used in the sample preservation, reaction of DNA amplification, electrophoresis, degeneration, solidification and other kinds of serum biochemical sample thermostatic incubation process.

Features:

- ◆ Real-time temperature display, the countdown display.
- ◆ Convenient module replacement, easy to clean and disinfect.
- ◆ Automatic fault detection and alarm function.
- ◆ Built-in over temperature protection device.
- ◆ Temperature deviation calibration function.

Parameters:

Project	Parameters
Temp.Range	RT.+5C -160
Temp uniformity	±0.5
Display accuracy	0.1
Heating time	≤15 min (25 to160)
Max power	500W
Size(mm)	260×220×100

Project	Parameters
Temp accuracy	±0.5 (@40)
Temp accuracy	±1 (@120)
Time range	1-99h59min/∞
Block quantity	2
Voltage	220V 50/60Hz
Weight(kg)	5.5

Optional Blocks:

Type	Diameter of test tube	Number of test tube
MD01	6mm	48
MD02	7mm	48
MD03	10mm	24
MD04	12mm	24
MD05	13mm	24
MD06	15mm	16
MD07	16mm	16
MD08	19mm	12
MD09	20mm	12

Type	Diameter of test tube	Number of test tube
MD10	26mm	8
MD11	28mm	4
MD12	40mm	3
MD13	0.2ml centrifuge tube	48
MD14	0.5ml centrifuge tube	48
MD15	1.5ml centrifuge tube	24
MD16	2.0ml centrifuge tube	24
MD17	0.2ml well micro plate	96
---	---	---

Novazym
life science & technology

Novazym Polska S.C. Wielkopolska Centre of Advanced Technologies
ul. Umultowska 89 C build. A/38, 61-614 Poznań