

Link do produktu: <http://www.novazym.sklep.pl/przeciwciala-poliklonalne-polyclonal-antibody-service-pakiet-i-p-439.html>



Przeciwciała poliklonalne (polyclonal antibody service) - PAKIET I

Dostępność

Dostępny

Numer katalogowy

PL001-01

Opis produktu

Custom Polyclonal Antibody Production:

This package is designed for generating antibodies against a client-supplied antigen. The package uses a client-provided antigen, and includes an 8-week immunization protocol in a species of your choice. We can perform services (ELISAs, purifications, western blotting, IHC etc.) for any project to assist you in evaluating the antibodies produced.

	Delivery	Total Time
Client-provided antigen	2-5mg lyophilized affinity-purified antibodies	15-16 weeks
2 rabbits or 1 goat	Lyophilized pre-immune serum	
Pre-Immune Serum from each animal	Elisa reports of 3 rd immunization serum	
4 immunizations per animal	Elisa reports of post-affinity antibodies	
Production Bleeds		
Peptide affinity purification		
Elisa test		