

Link do produktu: <http://www.novazym.sklep.pl/sample-concentrator-md200-1-p-377.html>



Sample Concentrator MD200-1

Dostępność

Na zamówienie

Numer katalogowy

MD200-1

Opis produktu

Sample Concentrator MD200-1

MD200-1 Sample Concentration consists of a base and stand assembly, sample holder and gas distribution system. Test tubes are held in a block .

Specification

1. Temperature range: RT. +5~150 °C
2. Temperature control discrepancy: $\leq \pm 0.5$ °C (RT +5 °C ~ 150 °C)
 $\leq \pm 1$ °C (100 °C ~ 150 °C)
3. Low Nitrogen consumption of 330 mL/min/sample
4. Time required for temperature to increase from 40 °C to 150 °C: ≤ 30 minutes
5. Maximum input power: 200 W
6. Temp. control discrepancy: ± 0.1 °C
7. Dimensions: 280W x 240D x 500H
8. Net weight: 6 kg
9. All units will fit and are recommended to be used in a bench-top.ventiated hood
10. Needles length: 150 mm
11. Standard Blocks: See "Blocks for MD200"

Key features

1. High precision and wide range of temperature control
2. Blowing a large number of samples at one time
3. Blowing each sample independently
4. Amount of airflow can be adjusted
5. Easy operation; safe and reliable to use

Blocks for MD200-1

Type	Parameter	Capacity	Dimension (mm)
BK01	6 mm	49	96 x 96 x 49
BK02	7 mm	49	96 x 96 x 49
BK03	10 mm	25	96 x 96 x 49
BK04	12 mm	25	96 x 96 x 49
BK05	13 mm	25	96 x 96 x 49
BK06	15 mm	12	96 x 96 x 49
BK07	15 mm	16	96 x 96 x 49
BK08	16 mm	12	96 x 96 x 49
BK09	16 mm	16	96 x 96 x 49
BK10	19 mm	12	96 x 96 x 49

BK11	19 mm	16	96 x 96 x 49
BK12	20 mm	9	96 x 96 x 49
BK13	26 mm	9	96 x 96 x 49
BK14	28 mm	4	96 x 96 x 49
BK15	40 mm	4	96 x 96 x 49
BK16	0.5 ml Centrifuge tube	49	96 x 96 x 49
BK17	1.5 ml Centrifuge tube	25	96 x 96 x 49
BK18	2.0 ml Centrifuge tube	25	96 x 96 x 49