

Link do produktu: <http://www.novazym.sklep.pl/termowyrzasarka-ms-100-thermo-shaker-incubator-p-133.html>



Termowyrząsarka MS-100 (Thermo Shaker Incubator)

Dostępność

Na zamówienie

Numer katalogowy

MS-100

Opis produktu

MS-100 Thermo-Shaker Incubator

Designed for simultaneous heating, cooling and mixing of small samples. Both MS-100 and MSC-100 are supplied with interchangeable platforms for various tubes.

Mixing, heating and cooling modes can be used either simultaneously or independently. The main body housing the mixing and heating platform can accommodate different kinds of blocks.

Features

- Clear, easy to use LCD display
- Soft touch keypads
- Small footprint
- Accurately control and display time, temperature and speed
- Over-heating protection device ensures safety & reliability
- Temperature can be calibrated to meet user's needs.
- Low noise working even under the speed of 1,500rpm
- Conforms to CE safety standards

Specification

Model	MS-100
Temp. Setting Range	5?~100?
Timing Range	1 min~99h59min
Mixing Speed	200rpm-1500rpm
Mixing Orbit	2mm
Temp. Control Accuracy	?? 0.3 ?
Display Accuracy	0.1 ?
Heating Time (from 20 ? to 100 ?)	?20min
Cooling time	-----
Dimension (LxWxH) mm	300x220x170

Net Weight

7 kg