

Link do produktu: <http://www.novazym.sklep.pl/id-card-centrifuge-td2-24-p-517.html>



ID Card Centrifuge TD2-24

Dostępność

Dostępny

Opis produktu

ID Card Centrifuge

Model: TD2-24

Digital centrifuge for Gel cards. The TD2-24 ID-Centrifuge holds automatic balance control and user can set speed and centrifugation time. It can be suitable many types ID-CARDS from different manufacturers.

Features

1. Simultaneous centrifuge of up to 24 Gel cards
2. Removable and exchangeable spinning head
3. Pre-defined centrifugal parameters for optimum Gel cards processing: acceleration, speed, braking and time
4. Digital interface visor control and visualization of centrifuge speed and time remaining for processing
5. Programmable audible alarms for end of centrifugal time periods

Specification

1. Capacity: 24 Gel cards
2. Speed: Two stages (It can be variable according to user need)
3. Max..Speed: 1500rpm
4. Timer: 0~99min
5. Power: AC100-240V, 50/60Hz, 50W
6. Dimensions: DxWxH=530x410x180 mm

Digital centrifuge for Gel cards. The TD2-24 ID-Centrifuge holds automatic balance control and user can set speed and centrifugation time. It can be suitable many types ID-CARDS from different manufacturers.

Features

1. Simultaneous centrifuge of up to 24 Gel cards
2. Removable and exchangeable spinning head
3. Pre-defined centrifugal parameters for optimum Gel cards processing: acceleration, speed, braking and time
4. Digital interface visor control and visualization of centrifuge speed and time remaining for processing

5. Programmable audible alarms for end of centrifugal time periods

Specification

1. Capacity: 24 Gel cards
2. Speed: Two stages (It can be vairable according to user need)
3. Max..Speed: 1500rpm
4. Timer: 0~99min
5. Power: AC100-240V, 50/60Hz, 50W
6. Demensions: DxWxH=530x410x180 mm