

Dane aktualne na dzień: 05-05-2024 05:17

Link do produktu: <http://www.novazym.sklep.pl/id-card-incubator-k37-24-p-516.html>



ID Card Incubator K37-24

Dostępność

Dostępny

Opis produktu

ID Card Incubator

Model: K37-24

Digital incubator for the Gel cards. The K37-24 ID-incubator holds 2 independent incubation zones with a capacity for 12 Gel cards, and user can set temperature and incubator time. It can be suitable many types ID-CARDS from different manufacturers.

Features

1. Incubator for 24-Gel Cards in racks at 37°C
2. 2 independent incubation zones with a capacity for 12 Gel cards
3. Digital interface screen for control and visualization of remaining incubation time and real temperature for each zone.
4. Fixed pre-set temperature of 37 degree
5. Pre-set incubation time of 15 minutes, which can be adjust by the user
6. Audible programmable alarm to notify end of incubation period

Specification

1. Capacity: 24 Gel cards
2. Temperature: 37°C
3. Timer: 15min (Variable)
4. Power: AC220V or AC120V 50-60Hz, 300W
5. Dimensions: DxWxH=380x220x115 mm