

Link do produktu: <http://www.novazym.sklep.pl/jy-zy6-western-blotting-system-trans-blot-core-p-362.html>

JY-ZY6 Western Blotting System (Trans-Blot Core)

Dostępność

Na zamówienie

Numer katalogowy

JY-ZY6

Opis produktu

JY-ZY6 Western Blotting System

Western blot system is used for testing a mixture of a great number of proteins. it can show us the information of our target protein (with comparison to a size marker or ladder in kDa), also the protein expression could be expressed by the run (with comparison to a control such as untreated sample or another cell type or tissue). "JY" series Western blot system covers almost all the areas of the general running, meeting the most people's experiments demands.

Advantages

JY-ZY6 Trans-Blot core is the main inner part of JY-ZY5, for the users who have purchased the JY-SCZ2 or JY-SCZ2+, they only need to equip the JY-ZY6 core to realize the transfer print electrophoresis run. JY-ZY6 Trans-Blot is equipped the identical with the JY-SCZ2+, customers can do the experiments just as the JY-SCZ2+ or JY-SCZ2.