

Dane aktualne na dzień: 04-06-2024 16:25

Link do produktu: <http://www.novazym.sklep.pl/jy02g-fast-gel-imaging-system-p-373.html>



JY02G Fast Gel Imaging System

Dostępność	Na zamówienie
Numer katalogowy	JY02G

Opis produktu

[JY02G Fast Gel Imaging System](#)

JY series UV Transilluminator is mainly applied to the electrophoresis results observation. With the principles of different wavelength shows the different electrophoresis band result, They are compatible with the observation of nuclear acid and protein electrophoresis bands, by using the high resolution of CCD and software technology, the gel pictures are clearly showed.

[Specification](#)

- Reflection ultraviolet wavelength: 254 nm 365 nm
- Transmission source ultraviolet wavelength: 312 nm
- Transmission size: 200 × 150 mm
- Features with high resolution of 1.4 million pixels and special optical filter
- Safety gel cutting door for operation, convenient and quick