

Dane aktualne na dzień: 14-05-2024 09:40

Link do produktu: <http://www.novazym.sklep.pl/jy5000-high-voltage-power-supply-p-342.html>



JY5000 High-voltage Power Supply

Dostępność	Na zamówienie
Numer katalogowy	JY5000

Opis produktu

[Specification](#)

- Size: 36×25×13(cm)
- Weight: 8.5kg
- Output Range: 20--5000V, 1--200mA, 1--100W.
- Output jacks: Two sets of output jacks are provided to facilitate connection of up to 2 identical electrophoresis cells simultaneously
- Timer control: 1min--99hr, 59min
- Display: Backlit LCD (128×64pixel)
- Control: Intellectual microprocessor control unit:

[Key features](#)

- Programmable constant voltage, constant current or constant
- Power operation with automatic crossover.
- Safety features include automatic detection of no-load, short circuit, rapid resistance change, ground leak and system overheating.
- Ten general electrophoresis methods stored function
- Automatic memory function
- Gel temperature testing function (need to purchase relevant temperature sensor separately)
- Standard timer operation function
- Automatic completion of a run interrupted by power failure
- Intellectual display of Constant-Voltage, Current and Power
- Optional functions (Recovery of power failure, Non-proliferation electric field, V-H?run step by step)