

Dane aktualne na dzień: 11-05-2024 10:51

Link do produktu: <http://www.novazym.sklep.pl/multi-speed-mini-centrifuge-hwlab-p-3806.html>



Multi speed mini centrifuge HWLAB

Opis produktu

Multi speed mini centrifuge HWLAB



Description

HWLAB® Multi speed mini centrifuge is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Product features

- Unique shape design
- LED speed display
- Adjustable speed from 1000r/min to 12000r/min
- Wide voltage, wide frequency

Technical specifications


Description	LED12k
Maximum rcf	6900xg
Speed	from 1000r/min to 12000r/min
Rotor capacity	6 x 1.5ml/2.0ml, 2 x 8 x 0.2ml, 4 x 5ml, 12 x 30mm/40mm
Adapter	6 x 0.2ml/0.5ml, 4 x 1.3ml
Power supply	AC100V-240V/50-60Hz
Dimension	17 x 15 x 13cm
Quantity per carton	4pcs
G.W.	6.0kg
Dimension per carton	44 x 40 x 36cm




Rotor #1



Rotor #2



Rotor #3



Rotor #4

Ordering information: HWLAB® Multi speed mini centrifuge

Model	Cat.NO	Destription
LED 12k	L-2101	Widely ues in clinical labs and donation centers
Rotor #1	L-2102	For 6x0.2ml/0.5ml/1.5ml/2.0ml tubes
Rotor #2	L-2103	For 2x8x0.2ml tubes
Rotor #3(Purchase)	L-2104	For 4x1.3ml/5ml tubes
Rotor #4(Purchase)	L-2105	For 12x30mm/40mm capillary tubes