



Link do produktu: <http://www.novazym.sklep.pl/clinical-centrifuge-hwlab-p-3812.html>



Opis produktu

Clinical centrifuge HWLAB





Description

HWLAB® Clinical centrifuge is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers.

Technical specifications

Description	HW12M
Maximum rcf	2460xg
Power supply	AC100V-240V/50-60Hz
Rotor capacity	12 x 15ml/10ml/7ml/5ml tubes or 8 x 50ml
Speed	From 300r/min to 4000r/min
Timer	15s-99min / continuous operation
Motor type	Brushless motor
Dimension	370x330x175mm
Weight	11kg

Product features

- Maintenance-free brushless motor
- LCD display of time and speed
- Adjustable speed from 300r/min to 4000r/min
- Timer, from 15s to 99min, or continuous operation
- For 12×10ml/7ml/5ml blood collection tubes or 12×15ml centrifuge tube or 8×50ml centrifuge tube
- Lid automatically open at the end of the run to prevent sample warming and to reduce handling time



Ordering information: HWLAB® Clinical centrifuge

	Model	Ord.No	Description
	HW12M	C-13000	Widely ues in clinical labs and donation centers
	HW12M.1	C-13001	For 12x15ml/10ml/7ml/5ml tubes
	HW12M.2	C-13002	For 8x50ml tubes