

Link do produktu: <http://www.novazym.sklep.pl/microplate-incubator-mk100-2a-p-306.html>



Microplate Incubator MK100-2A

Dostępność

Na zamówienie

Numer katalogowy

MK100-2A

Opis produktu

Microplate Incubator MK100-2A

Specification

The MK100-2A is a high performance microplate incubator which accommodates up to two microplates. The MK100-2A can be used for any enzyme or cell-based assays requiring uniform and strictly controlled incubation up to 70°C.

Incubation	Incubation range from 8 °C above ambient to 70 °C
	0.1°C steps
	Temperature stability ± 0.3 °C
	Temperature gradient < 0.5 °C across the plate (@ 37 °C)
Warming Up Speed	< 25 min (from 20 °C to 70 °C)
Capacity	Up to 2 microplates
	Max. microplate dimensions:
Power Requirements	115/230 V, 60/50 Hz, max. 150 VA
Dimensions	280 × 270 × 140 mm

Key features

- Microprocessor controlled time and temperature
- Beep-signal and stop of after program completion
- Heating plates above and below the microplates
- Excellent temperature accuracy and uniformity across the plate
- User definable parameters - temperature and time
- Platform capacity: 2 microtest plates