

Dane aktualne na dzień: 01-06-2024 22:01

Link do produktu: <http://www.novazym.sklep.pl/sample-concentrator-nd200-2-p-3799.html>



Sample Concentrator ND200-2

Opis produktu

Sample Concentrator

ND200-2



Description:

- ◆ ND200-2 sample concentrator is with microprocessor controlled and PID fuzzy controlled technique. Instead of traditional water bath device, the heat-conducting medium is high purity aluminum. The Test tubes are held in a block, and blowing a large number of samples at one time. Can be used in residue analysis, commodity inspection, food, environment, pharmacy and biological product and other industries.

Features:

- ◆ High precision and wide range of temp control.
- ◆ Blowing a large number of samples at one time.
- ◆ Blowing each sample independently.
- ◆ Amount of airflow can be adjusted.
- ◆ Easy operation; safe and reliable to use.
- ◆ 12-month warranty.

Parameters:

Project	Parameters
Temp. Range	RT.+5 -160
Temp. Uniformity	±0.5
Display accuracy	0.1
Heating time	≤15 min (25 TO 160)
Max gas usage	15L/min
Max gas pressure	0.1Mpa
Max power	500W
Size	260*220*475mm

Project	Parameters
Temp accuracy	±0.5 (@40)
Temp accuracy	±1 (@120)
Time Range	1-99h59min/∞
Max vertical travel	200mm
Needle length	150mm
Block qty	2
Voltage	220V 50/60Hz
Weight	7.5KG

Optional Blocks:

Type	Test tube diameter	Capacity
MD50	10mm	12
MD51	12mm	12
MD52	13mm	12
MD53	15mm	12

Type	Test tube diameter	Capacity
MD54	16mm	12
MD55	19mm	12
MD56	20mm	12
---	---	---

