

Link do produktu: <http://www.novazym.sklep.pl/termowyrzasarka-msc-100-thermo-shaker-incubator-p-134.html>



# Termowyrząsarka MSC-100 (Thermo Shaker Incubator)

Dostępność

Na zamówienie

Numer katalogowy

MSC-100

## Opis produktu

### MSC-100 Thermo-Shaker Incubator

Designed for simultaneous heating, cooling and mixing of small samples. Both MS-100 and MSC-100 are supplied with interchangeable platforms for various tubes. Mixing, heating and **cooling modes can be used either simultaneously or independently. The main body housing the mixing and heating platform can accommodate different kinds of blocks.**

#### Features:

- Clear, easy to use LCD display
- Soft touch keypads
- Small footprint
- Accurately control and display time, temperature and speed
- Over-heating protection device ensures safety & reliability
- Temperature can be calibrated to meet user's needs.
- Low noise working even under the speed of 1,500rpm
- Conforms to CE safty standards

#### Specification:

Model	MSC-100	
Temp. Setting Range	0?~100?	
Temp. Control Range	0?~100?	Ambient Temp.?20?
	4?~100?	Ambient Temp.?25?
	10?~100?	Ambient Temp.?30?
Timing Range	1 min~99h59min	
Mixing Speed	200rpm-1500rpm	
Mixing Orbit	2mm	
Temp. Control Accuracy	?? 0.3 ?	
Display Accuracy	0.1 ?	
Heating Time (from 20 ? to 100 ?)	?15min	
Cooling time	?30min (from RT. to RT.-20 ?)	
	?15min (from 100 ? to 20 ?)	
Dimension ( LxWxH) mm	300x220x170	

**Net Weight**

**8.5 kg**